

## Blue Mix Silicone Putty "Siligum" 200g (100g + 100g)



2 Part Formula for Mold Making

Rating: Not Rated Yet

### Price

Price:\$ 49.75

[Ask a question about this product](#)

## Blue Mix Silicone Putty "Siligum" - 2 Part Formula for Mold Making


Blue Mix Silicone - 2 Part Mold Making Putty  
Make your own detailed molds in minutes!

### Instructions

1. Simply mix equal parts of the blue and white agents.
2. Mix until the colour is uniform.
3. Press the original item that you want to make a mold into the silicone compound.
4. Let the silicone cure.
- 5 Gently remove the original.


You now have a mold that you can use dozens of times to make an exact replica of your original item!

Working time (including mixing time) = approximately 1 min 45 seconds  
Curing time = approximately 30 minutes




**Blue Mix**  
F-0117F (Silicone mold/ 25g + 25g)  
F-0117A (Silicone mold/100g +100g)  
F-0117M (Silicone mold/ 200g +200g)


Blue Mix silicone is so flexible that you can make a mold with dried piece and fired piece will be same size as original mold.




① Take same amount of white and blue agent.



② Mix them well.



③ While the mixed agent is soft, push into the mold. It will be cured after 30 minutes. Take it out carefully not to make the piece bend over.



④ Push Art Clay into the mold and dry it by a hair dryer for 30 mins, take it out from the mold.

Working time (include mixing time)	Approximately 1 min 45 seconds
Curing time (include mixing time)	Approximately 30 mins

(※Curing time differs slightly due to temperature and humidity.)

2 / 2