

Blue Mix Silicone Putty "Siligum" 200g (100g + 100g)



art lay Blue Mix Silicone
2 Part Molding Putty



art lay Blue Mix Silicone
2 Part Molding Putty

2 Part Formula for Mold Making

Rating: Not Rated Yet

Price

Price:\$ 49.75

[Ask a question about this product](#)

Blue Mix Silicone Putty "Siligum" - 2 Part Formula for Mold Making

Blue Mix Silicone - 2 Part Mold Making Putty
Make your own detailed molds in minutes!

Instructions

1. Simply mix equal parts of the blue and white agents.
2. Mix until the colour is uniform.
3. Press the original item that you want to make a mold into the silicone compound.
4. Let the silicone cure.
- 5 Gently remove the original.

You now have a mold that you can use dozens of times to make an exact replica of your original item!

Art Clay Products: Blue Mix Silicone Putty "Siligum" 200g (100g + 100g)

Working time (including mixing time) = approximately 1 min 45 seconds
Curing time = approximately 30 minutes

Blue Mix silicone is so flexible that you can make a mold with dried piece and fired piece will be same size as original mold.



① Take same amount of white and blue agent.



② Mix them well.



③ While the mixed agent is soft, push into the mold. It will be cured after 30 minutes. Take it out carefully not to make the piece bend over.

④ Push Art Clay into the mold and dry it by a hair dryer for 30 mins, take it out from the mold.

Blue Mix

- FC117F (Silicone mold/ 25g + 25g)
- FC117A (Silicone mold/100g +100g)
- FC117W (Silicone mold/ 200g +200g)

Working time (include mixing time)	Approximately 1 min 45 seconds
Curing time (include mixing time)	Approximately 30 mins

(※Curing time differs slightly due to temperature and humidity.)